

AJ5803 型电脑微量注射泵

AJ5803 MICRO-INFUSION SYRINGE PUMP



产品主要技术指标

- 工作环境温度: 5℃-40℃
- 工作电源交直流两用: AC 220V (AC110V)50Hz 或 DC 12V
- 功耗不大于: 25VA
- 可装载 20ml 50ml 100ml 一次性注射器
- 流量设定范围:
20ml 注射器 0.1ml/hr-150ml/hr
50ml 注射器 0.1ml/hr-300ml/hr
100ml 注射器 0.1ml/hr-300ml/hr
最小增量 0.1ml/hr
- 输液准确度: $\pm 5\%$ (含注射器精度)
- 有“残余量”、“阻塞”、“完成”以及“电池欠压”声光报警功能
- 机内电池自动应急切换, 完全充电后在 50ml/hr 流量可连续使用 4 小时以上
- 体积(L × B × H): 335 × 142 × 120mm
- 重量: 2.5kg
- Ambient Temperature: 5°C - 40°C
- Power Supply: AC 220V 50Hz or DC 12V (AC110V available if required)
- Power Consumption: Less than 25VA
- Applicable Syringes: 20, 50, 100ml disposable syringe
- Flow Rate Option:
20ml syringe: 0.1ml/hr - 150ml/hr in 0.1ml/h increments
50ml syringe: 0.1ml/hr - 300ml/hr in 0.1ml/h increments
100ml syringe: 0.1ml/hr - 300ml/hr in 0.1ml/h increments
- Infusion Accuracy: $\pm 5\%$ (including syringe)
- Audible and Visual Alarms of "NEAR EMPTY", "COMPLETE", "OCCLUSION", "LOW BATT"
- 4 hours continuously operation at 50ml/hr with built-in rechargeable Ni-MH battery
- Dimension(L × B × H): 335 × 142 × 120mm
- Weight: 2.5kg

AJ5806 型一次性使用输注泵

AJ5806 DISPOSABLE ANALGESIC PUMP

■ AJ5806 型一次性使用输注泵, 产品采用微机控制, 为机电一体化技术, 具有功能强, 操作简便, 并有声光报警及提示功能。该产品体积小, 便于携带, 为病人自控镇痛治疗提供了方便, 同时提高了药物治疗的实际效果。

■ AJ5806 Disposable Analgesic Pump adopts the technology of microcomputer control. The electromechanical unification technology enables the product to have many effective functions. The operation is simple and convenient; acousto-optic alarm and indicator are more patient friendly. Due to its small size, this product is easy to be carried by patient. This pump provides convenience for Patient-Controlled Analgesia (P.C.A.), and at the same time raises the actual effect of medicine treatment.

产品主要技术指标

- 产品电源采用 2 节 5 号电池
- 液盒容量为 100ml(260ml)
- 流量设定范围 1-10ml/hr
- 报警及提示:
a. 电池欠压提示功能 b. 报警提示(输液完成、运行故障)
c. 正常运行指示
- 有病人自控加注功能(PCA)
- 加注锁定时间不小于 30 分钟, 每次加注量为推进流量四分之一
- 体积(L × B × H): 140 × 95 × 35mm
- 重量: 230g
- Power supply: 2 AA batteries
- Medication Reservoir Capacity: 100ml or 260ml
- Flow Rate Options: Purge: 1, 2, 3, 4, 6, 8, 10ml/hr
- Audible and Visual Alarms of "COMPLETE / ERROR" and "LOW BATT"
- Patient-Controlled Analgesia (P.C.A.) Function: Each PCA dose is 1/4 of preset delivery dose. The interval between two P.C.A injections is no less than 30 minutes
- Dimensions (L × B × H): 140 × 95 × 35mm
- Weight: 230g

