


EOLO



MEFAR s.r.l. Via dei Prati, 62 - 25073 Bovezzo (BS) ITALY - Tel. 030/2713161 - 2003102 Fax 030/2000551

E-mail : medicare@satyam.net.in * www.medicareagencies.com

HIGH QUALITY AEROSOL MIST

EOLO aerosol therapy unit has been successfully used in the treatment of: BRONCHIAL ASTHMA - CHRONIC BRONCHITIS - EMPHYSEMA - MYCOSIS - BRONCHIECTASIS - CYSTIC FIBROSIS and other respiratory disorders.

Through the nebulizer MB2 or MB5, EOLO produces 2 cc of very fine mist in about 7 minutes with the particles in 0,5 to 5,5 micron range to ensure the most effective inhalation therapy for the entire respiratory tract.

MODERN STYLING AND PRACTICAL


EOLO has a compact and modern design with a storage compartment to hold accessories clean between the treatments.

It is simple to use with single-hand operations, ideal for paediatric and adult use at home and fully portable for the patients who want to travel without leaving their therapy behind.

RELIABLE

Each EOLO unit is manufactured under most exacting and controlled conditions to assure consistently high quality. Lubricant-free piston compressor allows quiet functioning and is virtually service-free.

A thermal over-load switch, available on request, protects the motor from overheating. Dust and impurities are kept by the foam filter placed on the body of the unit, as indicated on the drawing.



SPECIFICATIONS

Size mm: 165 L x 185 W x 275 H

Weight: 2.5 Kg


Voltage: on request

Units with particular performances, voltages and components are available on request.

ACCESSORIES

- PERMANENT NEBULIZER MB2 in three detachable parts: upper (1), lower half (2) and baffle (3), complete with "T" CONNECTOR (4) and MOUTH PIECE (5).
- NOSE PIECE (6).
- Non toxic 1 meter PVC MEDICAL TUBING (7), to connect the air outlet of the compressor to the lower half of nebulizer.
- ADULT MASK (8).
- PAEDIATRIC MASK (9) - optional.
- DISPOSABLE NEBULIZER MB5 - optional.

The set is complete and no extra accessories are needed for treatment.



WARRANTY

Each EOLO apparatus is covered with Mefar one year limited warranty if used in accordance with the instructions for use, on condition that all modifications or repairs are effected by persons authorized by Mefar and if connected to electrical installations which comply with the applicable national requirements.

Mefar shall therefore not be liable for any default, delay in performance caused by unauthorized tampering or by use not made according to instructions, by anomalies or faults of the electrical input line or by damages caused by rough handling during transport.

MANUFACTURED & MARKETING IN INDIA BY :



VENTILATOR
CPAP & BIPAP
PULSE OXIMETER
OXYGEN CONCENTRATOR
DEFIBRILLATOR WITH MONITOR

NEBULIZERS
CAPNOGRAPH
SLEEP LAB/SLEEP STUDY
PULSE OXIMETER & PROBES
MULTIPARAMETER MONITOR

203-205, VARDAN HOUSE, 7/28, ANSARI ROAD, DARYA GANJ,
NEW DELHI 110002 (INDIA)